

HARBORVIEW DRUG STICKER – *SAMPLES ONLY*

The pediatric emergency tape (BROSELOW™) shall be used to determine the approximate weight in kilograms and the corresponding color zone. Please see rules of use on the reverse side of these stickers.

Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL

Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL

Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL

Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL
Atropine 0.1 mg/mL	
PINK	1 mL
RED	1 mL
PURPLE	1 mL
YELLOW	1.3 mL
WHITE	1.6 mL
BLUE	2 mL
ORANGE	2.6 mL
GREEN	3.2 mL

Rules of Implementation

Implementation of this Color Coded Medication Sticker System requires 4 steps:

1. Double check medication concentration for accurate calculation of sticker dosages.
Different drug concentrations will give you different milliliter doses.
2. Discuss implementation with your pharmacy committee, code and pediatric committees, and emergency department supervisors to ensure agreement on specific doses.
3. Once dosages and sticker system has been approved, educate your staff during staff meetings, in-services days, skills days, etc.
4. Practice.
Especially in hospitals that rarely see children, we recommend quarterly pediatric mock codes where staff can practice using all pediatric resuscitation equipment, as well as practice medication dosages and delivery.